



		Femoral Sizes					
		C	D	E	F	G	H**
Tibial Size	1-2	CDEFGH/ 1-2 STRIPED PURPLE					
	3-4	CDEFGH/ 3-4 SOLID YELLOW					
	5-6	CDEFGH/ 5-6 SOLID GREEN					
	7-10	CDEFGH/ 7-10 SOLID BLUE					
Patella Size		Use standard size patellas with C-H standard femoral components			26mm (inset only)	32mm	38mm
					29mm	35mm	41mm

** CR only