



Zimmer® Periarticular Distal Radial Dorsal Locking Plates



Plate Features:

- Anatomically contoured plates create a fit that requires little or no additional bending and helps with metaphyseal/diaphyseal reduction
- The low profile plate facilitates fixation without impinging on soft tissues
- Plates are available in a variety of sizes and lengths, left and right
- Dual-compression slots will accommodate periarticular screws or conventional stainless steel screws and allow bi-directional compression
- The last diaphyseal plate hole is designed to accommodate the tension device (00-4817-000-05)
- The plate shaft design allows for a minimally invasive technique with submuscular passage of the plate
- Multiple locking holes in the plate head allow placement of the screws to capture fragments
- Thick-to-thin profiles make the plates autocontourable
- The anatomical shape of the head of the plate matches the shape of the dorsal surface of the distal radius

Prod. No. (Left)	Prod. No. (Right)	Description	Holes	Length
2.4mm Distal Dorsal Radial Locking Plates				
00-2358-022-02	00-2358-021-02	Distal Dorsal T Radial Locking Plate	2	36mm
00-2358-022-03	00-2358-021-03	Distal Dorsal T Radial Locking Plate	3	44mm
00-2358-022-04	00-2358-021-04	Distal Dorsal T Radial Locking Plate	4	56mm
00-2358-022-06	00-2358-021-06	Distal Dorsal T Radial Locking Plate	6	77mm
00-2358-022-08	00-2358-021-08	Distal Dorsal T Radial Locking Plate	8	98mm
00-2358-022-10	00-2358-021-10	Distal Dorsal T Radial Locking Plate	10	119mm
00-2358-022-12	00-2358-021-12	Distal Dorsal T Radial Locking Plate	12	139mm
00-2358-026-01	00-2358-025-01	Distal Dorsal Delta Radial Locking Plate	1	40mm
00-2358-026-02	00-2358-025-02	Distal Dorsal Delta Radial Locking Plate	2	52mm
00-2358-026-04	00-2358-025-04	Distal Dorsal Delta Radial Locking Plate	4	74mm
00-2358-026-06	00-2358-025-06	Distal Dorsal Delta Radial Locking Plate	6	95mm
00-2358-026-08	00-2358-025-08	Distal Dorsal Delta Radial Locking Plate	8	117mm
00-2358-026-10	00-2358-025-10	Distal Dorsal Delta Radial Locking Plate	10	138mm

Prod. No.	Description	Set Qty.
00-2360-000-23 Distal Radial Standard Jig Set		
00-2360-098-01	Dorsal T Jig, Right	1
00-2360-098-02	Dorsal T Jig, Left	1
00-2360-099-01	Dorsal Delta Jig, Right	1
00-2360-099-02	Dorsal Delta Jig, Left	1
00-2358-025-00	Distal Radial Plate and Jig Case	1

Prod. No.	Description	Set Qty.
00-2360-000-03 2.4mm Locking Instrument Set		
00-2358-045-00	Distal Radial Locking Screw and Instrument Case	1
00-2360-022-18	2.0mm/1.8mm Standard Cannula	4
00-2360-042-24	2.4mm/1.8mm Locking Screw Depth Gauge	1
00-2360-052-24	2.4mm Locking Screw Tap	2
00-2360-062-15	1.5mm Hex Forward Captive Screwdriver, QC	2
00-2360-062-25	2.5mm Hex Forward Captive Screwdriver, QC	2
00-2360-155-18	1.8mm Standard Drill	2
00-2360-175-20	2.0mm QC Drill Bit	2
00-4806-121-24	2.4mm QC Drill Bit (For gliding hole)	2
00-4808-024-01	2.4mm/1.8mm Double Drill Sleeve	1
00-4811-015-00	Screwdriver Handle, ZPS Mini, stainless steel,	1
00-4901-012-15	1.25mm Diameter K-wire, One end Trocar, 6" long	2
00-4901-016-15	1.6mm Diameter K-wire, One end Trocar, 6" long	2
00-2348-035-00	3.5mm Dual Compression Drill Guide	1
00-4806-110-25	2.5mm Drill Bit, 110mm	2
00-4806-110-35	3.5mm Drill Bit, 110mm (For gliding hole)	2
00-4808-035-01	3.5mm/2.5mm Double Drill Sleeve	1
00-4808-035-02	3.5mm/2.5mm Drill Sleeve	1
00-4810-002-01	Small Depth Gauge	1
00-4811-035-00	T-Handle, QC	1
00-4811-110-35	3.5mm Tap	2
00-4812-001-00	Screw Forceps	1
00-4812-035-00	Small Hex Screwdriver, 2.5mm Hex	1
00-4812-035-05	Small Holding Sleeve	1
00-4819-051-00	Pin Cutter (Used to cut 1.8mm Locking Pegs)	1
00-4812-015-05	Mini Holding Sleeve	1
00-4812-015-01	Mini Hex Screwdriver Shaft	1

Designed for use with:

Zimmer 2.4mm Periarticular Locking Screw Set (00-2359-000-03)
Zimmer 2.4mm Periarticular Locking Instrument Set (00-2360-000-03)
Zimmer Periarticular Distal Radial Locking Plate Standard Jig Set (00-2360-000-23)

Compatible with:

Zimmer Periarticular 3.5mm Self-tapping Cortical Screws with 2.7mm Head
Zimmer 2.7mm Mini Fragment Cortical Screws
Zimmer Small Fragment Instrument Set (00-4800-035-20)





Prod. No.	Description	Set Qty.
00-2359-000-03	2.4mm Locking Screw Set	
00-2359-008-24	2.4mm Locking Screw 8mm Lng	4
00-2359-010-24	2.4mm Locking Screw 10mm Lng	4
00-2359-012-24	2.4mm Locking Screw 12mm Lng	4
00-2359-014-24	2.4mm Locking Screw 14mm Lng	6
00-2359-016-24	2.4mm Locking Screw 16mm Lng	6
00-2359-018-24	2.4mm Locking Screw 18mm Lng	6
00-2359-020-24	2.4mm Locking Screw 20mm Lng	6
00-2359-022-24	2.4mm Locking Screw 22mm Lng	4
00-2359-024-24	2.4mm Locking Screw 24mm Lng	4
00-2359-026-24	2.4mm Locking Screw 26mm Lng	4
00-2359-028-24	2.4mm Locking Screw 28mm Lng	4
00-2359-030-24	2.4mm Locking Screw 30mm Lng	4
00-2359-032-24	2.4mm Locking Screw 32mm Lng	2
00-2359-034-24	2.4mm Locking Screw 34mm Lng	2
00-2359-036-24	2.4mm Locking Screw 36mm Lng	2
00-2359-038-24	2.4mm Locking Screw 38mm Lng	2
00-2359-040-24	2.4mm Locking Screw 40mm Lng	2
00-2359-018-18	1.8mm Locking Peg 18mm Lng	4
00-2359-020-18	1.8mm Locking Peg 20mm Lng	4
00-2359-022-18	1.8mm Locking Peg 22mm Lng	4
00-2359-024-18	1.8mm Locking Peg 24mm Lng	4
00-2359-026-18	1.8mm Locking Peg 26mm Lng	4
00-2359-028-18	1.8mm Locking Peg 28mm Lng	4
00-2359-030-18	1.8mm Locking Peg 30mm Lng	4
00-2359-040-18	1.8mm Locking Peg 40mm Lng	4
00-2359-008-25	2.4mm Conical Screw 8mm Lng	4
00-2359-010-25	2.4mm Conical Screw 10mm Lng	4
00-2359-012-25	2.4mm Conical Screw 12mm Lng	4
00-2359-014-25	2.4mm Conical Screw 14mm Lng	6
00-2359-016-25	2.4mm Conical Screw 16mm Lng	6
00-2359-018-25	2.4mm Conical Screw 18mm Lng	6
00-2359-020-25	2.4mm Conical Screw 20mm Lng	6
00-2359-022-25	2.4mm Conical Screw 22mm Lng	4
00-2359-024-25	2.4mm Conical Screw 24mm Lng	4
00-2359-026-25	2.4mm Conical Screw 26mm Lng	4
00-2359-028-25	2.4mm Conical Screw 28mm Lng	4
00-2359-030-25	2.4mm Conical Screw 30mm Lng	4
00-2359-032-25	2.4mm Conical Screw 32mm Lng	2
00-2359-034-25	2.4mm Conical Screw 34mm Lng	2
00-2359-036-25	2.4mm Conical Screw 36mm Lng	2
00-2359-038-25	2.4mm Conical Screw 38mm Lng	2
00-2359-040-25	2.4mm Conical Screw 40mm Lng	2



Prod. No.	Description	Set Qty.
Existing Screws Included In Set No. 00-2359-000-03		
00-2348-008-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 8mm Lng	4
00-2348-010-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 10mm Lng	4
00-2348-012-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 12mm Lng	6
00-2348-014-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 14mm Lng	6
00-2348-016-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 16mm Lng	6
00-2348-018-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 18mm Lng	6
00-2348-020-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 20mm Lng	6
00-2348-022-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 22mm Lng	6
00-2348-024-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 24mm Lng	6
00-2348-026-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 26mm Lng	6
00-2348-028-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 28mm Lng	6
00-2348-030-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 30mm Lng	6
00-2348-032-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 32mm Lng	4
00-2348-034-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 34mm Lng	4
00-2348-036-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 36mm Lng	4
00-2348-038-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 38mm Lng	4
00-2348-040-35	3.5mm Cortical Self-Tapping (2.7mm Head) Screw 40mm Lng	4



97-2347-106-00 10MM Printed in USA ©2005 Zimmer, Inc.

Contact your Zimmer representative or
visit us at www.zimmer.com

