



Zimmer® Periarticular Proximal Lateral Humeral Locking Plate

Plate Features:

- Anatomically contoured plates create a fit that requires little or no additional bending and are helpful for metaphyseal/diaphyseal reduction
- The low profile plate facilitates fixation without impinging on soft tissues
- Plates are available in a variety of sizes and lengths, from 4 to 16 holes, left and right
- Dual-compression slots accommodate periarticular screws or conventional stainless steel screws and allow bi-directional compression
- The last diaphyseal plate hole is designed to accommodate the tension device (00-4817-000-05)
- The tapered plate shaft design allows for a minimally invasive technique with submuscular passage of the plate
- Multiple locking holes in the plate head allow placement of the screws to capture fragments
- Thick-to-thin profiles make the plates autocontourable
- The anatomical shape of the head of the plate matches the shape of the proximal humerus

Prod. No. (Left)	Prod. No. (Right)	Description	Holes	Length
3.5mm Proximal Lateral Humeral Locking Plate				
00-2358-002-04	00-2358-001-04	Proximal Lateral Humeral Locking Plate	4	90mm
00-2358-002-06	00-2358-001-06	Proximal Lateral Humeral Locking Plate	6	114mm
00-2358-002-08	00-2358-001-08	Proximal Lateral Humeral Locking Plate	8	138mm
00-2358-002-12	00-2358-001-12	Proximal Lateral Humeral Locking Plate	12	186mm
00-2358-002-16	00-2358-001-16	Proximal Lateral Humeral Locking Plate	16	234mm

Prod. No.	Description	Set Qty.
00-2358-000-05 Proximal Lateral Humeral Locking Plate Set		
00-2358-002-04	Proximal Lateral Humeral Locking Plate, 4 Hole 90mm Lng, Left	1
00-2358-002-06	Proximal Lateral Humeral Locking Plate, 6 Hole 114mm Lng, Left	1
00-2358-002-08	Proximal Lateral Humeral Locking Plate, 8 Hole 138mm Lng, Left	1
00-2358-002-12	Proximal Lateral Humeral Locking Plate, 12 Hole 186mm Lng, Left	1
00-2358-002-16	Proximal Lateral Humeral Locking Plate, 16 Hole 234mm Lng, Left	1
00-2358-001-04	Proximal Lateral Humeral Locking Plate, 4 Hole 90mm Lng, Right	1
00-2358-001-06	Proximal Lateral Humeral Locking Plate, 6 Hole 114mm Lng, Right	1
00-2358-001-08	Proximal Lateral Humeral Locking Plate, 8 Hole 138mm Lng, Right	1
00-2358-001-12	Proximal Lateral Humeral Locking Plate, 12 Hole 186mm Lng, Right	1
00-2358-001-16	Proximal Lateral Humeral Locking Plate, 16 Hole 234mm Lng, Right	1

Prod. No.	Description	Set Qty.
00-2360-000-14 Proximal Lateral Humeral Plate Standard Jig Set		
00-2360-092-01	Proximal Lateral Humeral Plate Jig, Right	1
00-2360-092-02	Proximal Lateral Humeral Plate Jig, Left	1
00-2358-030-00	Case	1

Prod. No.	Description	Set Qty.
00-2360-000-02 3.5mm/2.7mm Periarticular Locking Instrument Set		
00-1147-071-00	Cleaning Stylet	1
00-1147-076-00	Cleaning Brush	1
00-2358-040-00	3.5/2.7mm Locking Screw and Instrument Case	1
00-2360-011-01	3.5mm/2.7mm Locking Plate Reduction Instrument	2
00-2360-011-02	3.5mm/2.7mm Locking Plate Reduction Sleeve	2
00-2360-012-03	Plate Reduction Spin Knob	2
00-2360-020-20	2.0mm Standard Cannula	4
00-2360-020-27	2.7mm Standard Cannula	4
00-2360-021-16	1.6mm Standard Cannula	4
00-2360-033-16	1.6mm Standard Drill Tip Guide Wire (5 per box)	2 bx
00-2360-040-35	3.5mm/2.7mm Locking Screw Standard Depth Gauge	1
00-2360-041-35	3.5mm Cannulated Locking Screw Depth Gauge	1
00-2360-053-27	2.7mm Locking Screw Tap	2
00-2360-053-35	3.5mm Locking Screw Tap	2
00-2360-054-35	3.5mm Cannulated Locking Screw Tap	2
00-2360-065-00	2.5mm Screwdriver Stop Ring	1
00-2360-065-25	2.5mm Hex Standard Screwdriver	2
00-2360-066-25	2.5mm Hex Standard Cannulated Screwdriver	2
00-2360-071-27	2.7mm Standard Cannulated Drill	2
00-2360-080-00	Torque Limiting Attachment	1
00-2360-086-00	Modular Handle	1
00-2360-088-00	3.5mm/2.7mm Cannula Insert	1
00-2360-175-20	2.0mm Standard Drill	2
00-2360-205-27	2.7mm Standard Drill	2
00-4812-045-00	Large Hex Screwdriver	1
Also Available:		
00-2360-050-25	2.5mm Hex Forward Captive Screwdriver	0
00-2360-051-25	2.5mm Hex Cannulated Forward Captive Screwdriver	0

Designed for use with:

The 3.5mm screws in the *Zimmer* 3.5mm/2.7mm Periarticular Locking Screw Set (00-2359-000-02)
 The 3.5mm screws in the *Zimmer* 3.5mm/2.7mm Periarticular Locking Instrument Set (00-2360-000-02)
Zimmer Periarticular Proximal Lateral Humeral Locking Plate Standard Jig Set (00-2360-000-14)

Compatible with:

Zimmer Periarticular 3.5mm Self-tapping Cortical Screws
Zimmer Periarticular 4.0mm Cancellous Screws
Zimmer Small Fragment Instrument Set (00-4800-035-20)
Zimmer Small Fragment Screw Set (00-4800-035-01)





Prod. No.	Description	Set Qty.		
00-2359-000-02	3.5mm/2.7mm Locking Screw Set			
00-2359-030-36	3.5mm Cannulated Locking Screw 30mm Lng	3	00-2359-040-35	3.5mm Locking Screw 40mm Lng 3
00-2359-035-36	3.5mm Cannulated Locking Screw 35mm Lng	3	00-2359-042-35	3.5mm Locking Screw 42mm Lng 3
00-2359-040-36	3.5mm Cannulated Locking Screw 40mm Lng	3	00-2359-044-35	3.5mm Locking Screw 44mm Lng 3
00-2359-045-36	3.5mm Cannulated Locking Screw 45mm Lng	3	00-2359-046-35	3.5mm Locking Screw 46mm Lng 3
00-2359-050-36	3.5mm Cannulated Locking Screw 50mm Lng	4	00-2359-048-35	3.5mm Locking Screw 48mm Lng 3
00-2359-055-36	3.5mm Cannulated Locking Screw 55mm Lng	4	00-2359-050-35	3.5mm Locking Screw 50mm Lng 3
00-2359-060-36	3.5mm Cannulated Locking Screw 60mm Lng	4	00-2359-052-35	3.5mm Locking Screw 52mm Lng 2
00-2359-065-36	3.5mm Cannulated Locking Screw 65mm Lng	4	00-2359-054-35	3.5mm Locking Screw 54mm Lng 2
00-2359-070-36	3.5mm Cannulated Locking Screw 70mm Lng	4	00-2359-056-35	3.5mm Locking Screw 56mm Lng 2
00-2359-075-36	3.5mm Cannulated Locking Screw 75mm Lng	3	00-2359-058-35	3.5mm Locking Screw 58mm Lng 2
00-2359-080-36	3.5mm Cannulated Locking Screw 80mm Lng	3	00-2359-060-35	3.5mm Locking Screw 60mm Lng 2
00-2359-085-36	3.5mm Cannulated Locking Screw 85mm Lng	3	00-2359-065-35	3.5mm Locking Screw 65mm Lng 2
00-2359-090-36	3.5mm Cannulated Locking Screw 90mm Lng	3	00-2359-070-35	3.5mm Locking Screw 70mm Lng 2
00-2359-030-37	3.5mm Cannulated Conical Screw 30mm Lng	2	00-2359-010-27*	2.7mm Locking Screw 10mm Lng 3
00-2359-035-37	3.5mm Cannulated Conical Screw 35mm Lng	2	00-2359-012-27*	2.7mm Locking Screw 12mm Lng 3
00-2359-040-37	3.5mm Cannulated Conical Screw 40mm Lng	2	00-2359-014-27*	2.7mm Locking Screw 14mm Lng 3
00-2359-045-37	3.5mm Cannulated Conical Screw 45mm Lng	2	00-2359-016-27*	2.7mm Locking Screw 16mm Lng 3
00-2359-050-37	3.5mm Cannulated Conical Screw 50mm Lng	2	00-2359-018-27*	2.7mm Locking Screw 18mm Lng 3
00-2359-055-37	3.5mm Cannulated Conical Screw 55mm Lng	2	00-2359-020-27*	2.7mm Locking Screw 20mm Lng 3
00-2359-060-37	3.5mm Cannulated Conical Screw 60mm Lng	2	00-2359-022-27*	2.7mm Locking Screw 22mm Lng 4
00-2359-065-37	3.5mm Cannulated Conical Screw 65mm Lng	2	00-2359-024-27*	2.7mm Locking Screw 24mm Lng 4
00-2359-070-37	3.5mm Cannulated Conical Screw 70mm Lng	2	00-2359-026-27*	2.7mm Locking Screw 26mm Lng 4
00-2359-012-35	3.5mm Locking Screw 12mm Lng	2	00-2359-028-27*	2.7mm Locking Screw 28mm Lng 4
00-2359-014-35	3.5mm Locking Screw 14mm Lng	2	00-2359-030-27*	2.7mm Locking Screw 30mm Lng 4
00-2359-016-35	3.5mm Locking Screw 16mm Lng	2	00-2359-032-27*	2.7mm Locking Screw 32mm Lng 4
00-2359-018-35	3.5mm Locking Screw 18mm Lng	2	00-2359-034-27*	2.7mm Locking Screw 34mm Lng 4
00-2359-020-35	3.5mm Locking Screw 20mm Lng	2	00-2359-036-27*	2.7mm Locking Screw 36mm Lng 4
00-2359-022-35	3.5mm Locking Screw 22mm Lng	2	00-2359-038-27*	2.7mm Locking Screw 38mm Lng 4
00-2359-024-35	3.5mm Locking Screw 24mm Lng	2	00-2359-040-27*	2.7mm Locking Screw 40mm Lng 4
00-2359-026-35	3.5mm Locking Screw 26mm Lng	2	00-2359-042-27*	2.7mm Locking Screw 42mm Lng 4
00-2359-028-35	3.5mm Locking Screw 28mm Lng	2	00-2359-044-27*	2.7mm Locking Screw 44mm Lng 4
00-2359-030-35	3.5mm Locking Screw 30mm Lng	2	00-2359-046-27*	2.7mm Locking Screw 46mm Lng 4
00-2359-032-35	3.5mm Locking Screw 32mm Lng	3	00-2359-048-27*	2.7mm Locking Screw 48mm Lng 4
00-2359-034-35	3.5mm Locking Screw 34mm Lng	3	00-2359-050-27*	2.7mm Locking Screw 50mm Lng 4
00-2359-036-35	3.5mm Locking Screw 36mm Lng	3	00-2359-052-27*	2.7mm Locking Screw 52mm Lng 3
00-2359-038-35	3.5mm Locking Screw 38mm Lng	3	00-2359-054-27*	2.7mm Locking Screw 54mm Lng 3
			00-2359-056-27*	2.7mm Locking Screw 56mm Lng 3
			00-2359-058-27*	2.7mm Locking Screw 58mm Lng 3
			00-2359-060-27*	2.7mm Locking Screw 60mm Lng 3

* The proximal lateral humeral locking plate is designed for use with 3.5mm screws, and is not intended for use with the 2.7mm screw.



97-2347-108-00 10MM Printed in USA ©2005 Zimmer, Inc.

Contact your Zimmer representative or
visit us at www.zimmer.com

