



Cross Section Taken at Y-Y and C-C as indicated below



Cross Section Taken at Y-Y and C-C as indicated below



Catalog Number	Component Description	Length (mm) (in)	Cross Section Y-Y		Cross Section C-C	
			A (mm) x B (mm)	A (in) x B(in)	A (mm) x B (mm)	A (in) x B (in)
EXTRA-SMALL ELBOW						
32-8105-93-01	Left Ulnar Component	114.3 / 4.5	3.5 / 3.5	0.14 / 0.14	3.5 / 3.5	0.14 / 0.14
32-8105-93-02	Right Ulnar Component	114.3 / 4.5	3.5 / 3.5	0.14 / 0.14	3.5 / 3.5	0.14 / 0.14
32-8105-27-04	Humeral Component	101.6 / 4.0	6.6 / 4.5	0.26 / 0.18	7.1 / 4.5	0.28 / 0.18
32-8105-27-06	Humeral Component	152.4 / 6.0	5.6 / 4.5	0.22 / 0.18	7.1 / 4.5	0.28 / 0.18

SMALL ELBOW

32-8105-53-01	Left Ulnar Component	76.2 / 3.0	5.1 / 5.6	0.20 / 0.22	5.6 / 5.6	0.22 / 0.22
32-8105-53-02	Right Ulnar Component	76.2 / 3.0	5.1 / 5.6	0.20 / 0.22	5.6 / 5.6	0.22 / 0.22
32-8105-73-01	Left Ulnar Component	114.3 / 4.5	5.1 / 5.6	0.20 / 0.22	5.6 / 5.6	0.22 / 0.22
32-8105-73-02	Right Ulnar Component	114.3 / 4.5	5.1 / 5.6	0.20 / 0.22	5.6 / 5.6	0.22 / 0.22
32-8105-25-04	Humeral Component	101.6 / 4.0	7.3 / 5.1	0.29 / 0.20	7.9 / 5.1	0.31 / 0.20
32-8105-25-06	Humeral Component	152.4 / 6.0	6.4 / 5.1	0.25 / 0.20	7.9 / 5.1	0.31 / 0.20
32-8105-25-08	Humeral Component	203.2 / 8.0	5.6 / 5.1	0.22 / 0.20	7.9 / 5.1	0.31 / 0.20

REGULAR ELBOW

32-8105-63-01	Left Ulnar Component	88.9 / 3.5	5.8 / 5.8	0.23 / 0.23	6.4 / 5.8	0.25 / 0.23
32-8105-63-02	Right Ulnar Component	88.9 / 3.5	5.8 / 5.8	0.23 / 0.23	6.4 / 5.8	0.25 / 0.23
32-8105-83-01	Left Ulnar Component	114.3 / 4.5	5.6 / 5.8	0.22 / 0.23	6.4 / 5.8	0.25 / 0.23
32-8105-83-02	Right Ulnar Component	114.3 / 4.5	5.6 / 5.8	0.22 / 0.23	6.4 / 5.8	0.25 / 0.23
32-8105-26-04	Humeral Component	101.6 / 4.0	7.1 / 6.1	0.37 / 0.24	9.9 / 6.1	0.39 / 0.24
32-8105-26-06	Humeral Component	152.4 / 6.0	8.1 / 6.1	0.32 / 0.24	9.9 / 6.1	0.39 / 0.24
32-8105-26-08	Humeral Component	203.2 / 8.0	6.6 / 6.1	0.26 / 0.24	9.9 / 6.1	0.39 / 0.24

POCKET REFERENCE GUIDE

