

Reference Combination Manual/Automated Cleaning and Disinfection Instructions (Section H) of the Instructions for Care, Cleaning, Maintenance and Sterilization Manual (97-5000-170-00). Do not utilize the optional ultrasonic cleaning note in Section H.



Persona® Shims

42-5279-000-01_04

42-5279-002-01_04


42-5279-004-01_04

42-5279-006-01_04

42-5279-008-01-04

Cleaning Instructions

Do not ultrasonic clean



97-5026-055-00 Rev. 1 Printed in USA ©2015 Zimmer, Inc.